CE IVD  
REF 2041-3040

## SARS-CoV-2 Antigen Test Kit(Colloidal Gold)

## Instructions for Use

## INTENDED USE

The SARS-CoV-2 Antigen Test Kit is a gold immuno-chromatographic assay (GICA) that is intended for the qualitative detection of the nucleocapsid protein antigen from SARS-CoV-2 in nasopharyngeal (NP) swab, nasal (NS) swab, oropharyngeal (OP) swab and saliva specimens directly from individuals who are suspected of COVID-19 by their healthcare provider.

The SARS-CoV-2 Antigen Test Kit does not differentiate between SARS-CoV and SARS-CoV-2.

Results are for the detection of SARS-CoV-2 nucleocapsid protein antigen. Antigen is generally detectable in upper respiratory specimens during the acute phase of infection. Positive results indicate the presence of viral antigens, but clinical correlation with patient history and other diagnostic information is necessary to determine infection status. Positive results do not rule out bacterial co-infection or co-infection with other viruses.

The agent detected may not be definitive cause of disease.

Negative results should be treated as presumptive and confirmed with a molecular assay, if necessary for patient management. Negative results do not rule out COVID-19 and should not be used as the sole basis for treatment and patient management decisions, including infection control decisions. Negative results should be considered in the context of a patient's recent exposures, history and the presence of clinical signs and symptoms consistent with COVID-19.

Individuals who test negative and continue to experience COVID-like symptoms should seek follow up care from their healthcare provider.

## GENERAL INFORMATION

The novel coronaviruses belong to the *β* genus. COVID-19 is an acute respiratory infectious disease. People are generally susceptible. Currently, the patients infected by the novel coronavirus are the main source of infection; asymptomatic infected people also can be an infectious source. Based on the current epidemiological investigation, the incubation period is 1 to 14 days, mostly 3 to 7 days. The main manifestations include fever, fatigue and dry cough. Nasal congestion, runny nose, sore throat, myalgia and diarrhea are found in a few cases.

The median incubation time is estimated to be 5.1 days with symptoms expected to be present within 12 days of infection. The symptoms of COVID-19 are similar to other viral respiratory diseases and include fever, cough and shortness of breath.

## PRINCIPLE OF THE TEST

The SARS-CoV-2 Antigen Test Kit is a rapid lateral flow immuno-chromatographic sandwich assay to directly detect nucleocapsid protein of SARS-CoV-2 in nasopharyngeal swab, nasal swab, oropharyngeal (OP) swab and saliva specimens and diagnosis of SARS-CoV-2 infection.

The patient sample is placed in the Sample Tube, during which time the virus particles in the sample are disrupted, exposing internal viral nucleoproteins. After disruption, the sample is added into the test cassette sample well. And the sample migrates through a test strip. If the SARS-CoV-2 viral antigen is present, a red color line will be shown on the T line. If no SARS-CoV-2 viral antigen is absent, there is not a red line will be shown on the T line, however, a red line will always be showed on the C line indicating that the reaction system is properly happened.

## REAGENTS AND MATERIALS PROVIDED

Item	Component	Specifications/Qty.
1	Test Cassette individually foil packed with a desiccant	REF 203-001 / REF 203-020
2	SampleTube, with,0.5mLsample buffer	1 20
3	Single packaged swab(NS/OP/Pswab)*	1 20
4	Instructions for use	1 1
5*	Saliva collection device/bag with a 0.5mLTransfer Pipette*	1 20

\* Components will be included on customer's demand.

Materials needed but not provided:

1. Timer or watch,Vortex.

2. Please read this manual carefully prior to using this test kit. And follow the testing procedures strictly described in the manual, otherwise it will lead to incorrect results.

3. Do not use expired reagents.

4. Do not reuse the test kit.

5. All swab samples, tested reagents, tested cassettes and other materials used during testing are considered to be infectious, and personal protection should be done during the experiment.

6. Use of Nitrile Latex(equivalent) gloves is recommended when handling patient samples. Wear suitable protective clothing and eye/face protection when handling the contents of this kit.

7. Sample handling and waste disposal must comply with relevant regulations. Wash hands thoroughly after handling.

8. Avoid using visually bloody or overly viscous samples for testing.

9. Do not use components from different batch lots.

10. The sample tube contains a salt solution. If the solution contacts the skin or eye, flush with copious amounts of water.

11. Sample collection and handling procedures require specific training and guidance

## STORAGE AND STABILITY

1. The test device is sensitive to humidity as well as to heat.

2. Store kit components at 2-30°C, out of direct sunlight. Kit components are stable until the expiration date printed on the outer box.

shown on the T line, however, a red line will always be showed on the C line indicating that the reaction system is properly happened.

## REAGENTS AND MATERIALS PROVIDED

Item	Component	Specifications/Qty.
1	Test Cassette individually foil packed with a desiccant	REF 203-001 / REF 203-020
2	SampleTube, with,0.5mLsample buffer	1 20
3	Single packaged swab(NS/OP/Pswab)*	1 20
4	Instructions for use	1 1
5*	Saliva collection device/bag with a 0.5mLTransfer Pipette*	1 20

\* Components will be included on customer's demand.

Materials needed but not provided:

1. Timer or watch,Vortex.

2. Please read this manual carefully prior to using this test kit. And follow the testing procedures strictly

described in the manual, otherwise it will lead to incorrect results.

3. Do not use expired reagents.

4. Do not reuse the test kit.

5. All swab samples, tested reagents, tested cassettes and other materials used during testing are considered to be infectious, and personal protection should be done during the experiment.

6. Use of Nitrile Latex(equivalent) gloves is recommended when handling patient samples. Wear suitable protective clothing and eye/face protection when handling the contents of this kit.

7. Sample handling and waste disposal must comply with relevant regulations. Wash hands thoroughly after handling.

8. Avoid using visually bloody or overly viscous samples for testing.

9. Do not use components from different batch lots.

10. The sample tube contains a salt solution. If the solution contacts the skin or eye, flush with copious amounts of water.

11. Sample collection and handling procedures require specific training and guidance

## TESTPRINZIPI

The SARS-CoV-2 Antigen Test Kit is a CE marked RT-PCR test kit were evaluated.

SARS-CoV-2 Antigen Test Performance against the RT-PCR Method with Ascorpharyngeal/ Oropharyngeal swab Specimen

SARS-CoV-2 Antigen Test Kit	Comparator Method (RT-PCR)	Positive	Negative	Total
		114	2	116
		6	218	224
		120	220	340

Sensitivity(95% CI): 14/120 (95.0%-98.1%)

Specificity(95% CI): 218/220 (96.2%-99.5%)

Accuracy(95% CI): 114/218 (340) (97.5% (95.0%-99.0%))

In the clinical study of the SARS-CoV-2 Antigen Test Kit on nasopharyngeal/oropharyngeal swabs, 120 positive specimens and 220 negative specimens confirmed by a CE marked RT-PCR test kit were evaluated.

SARS-CoV-2 Antigen Test Performance against the RT-PCR Method with Ascorpharyngeal/ Oropharyngeal swab Specimen

SARS-CoV-2 Antigen Test Kit	Comparator Method (RT-PCR)	Positive	Negative	Total
		114	2	116
		6	218	224
		120	220	340

Sensitivity(95% CI): 14/120 (95.0%-98.1%)

Specificity(95% CI): 218/220 (96.2%-99.5%)

Accuracy(95% CI): 114/218 (340) (97.5% (95.0%-99.0%))

In the clinical study of the SARS-CoV-2 Antigen Test Kit on nasal swabs, 110 positive specimens and 220 negative specimens confirmed by a CE marked RT-PCR test kit were evaluated.

SARS-CoV-2 Antigen Test Performance against the RT-PCR Method with Nasal Swab Specimen

SARS-CoV-2 Antigen Test Kit	Comparator Method (RT-PCR)	Positive	Negative	Total
		114	2	116
		6	218	224
		120	220	340

Sensitivity(95% CI): 14/120 (95.0%-98.1%)

Specificity(95% CI): 218/220 (96.2%-99.5%)

Accuracy(95% CI): 114/218 (340) (97.5% (95.0%-99.0%))

In the clinical study of the SARS-CoV-2 Antigen Test Kit on nasal swabs, 110 positive specimens and 220 negative specimens confirmed by a CE marked RT-PCR test kit were evaluated.

SARS-CoV-2 Antigen Test Performance against the RT-PCR Method with Nasal Swab Specimen

SARS-CoV-2 Antigen Test Kit	Comparator Method (RT-PCR)	Positive	Negative	Total
		114	2	116
		6	218	224
		120	220	340

Sensitivity(95% CI): 14/120 (95.0%-98.1%)

Specificity(95% CI): 218/220 (96.2%-99.5%)

Accuracy(95% CI): 114/218 (340) (97.5% (95.0%-99.0%))

In the clinical study of the SARS-CoV-2 Antigen Test Kit on nasal swabs, 110 positive specimens and 220 negative specimens confirmed by a CE marked RT-PCR test kit were evaluated.

SARS-CoV-2 Antigen Test Performance against the RT-PCR Method with Nasal Swab Specimen

SARS-CoV-2 Antigen Test Kit	Comparator Method (RT-PCR)	Positive	Negative	Total
		114	2	116
		6	218	224
		120	220	340

Sensitivity(95% CI): 14/120 (95.0%-98.1%)

Specificity(95% CI): 218/220 (96.2%-99.5%)

Accuracy(95% CI): 114/218 (340) (97.5% (95.0%-99.0%))

In the clinical study of the SARS-CoV-2 Antigen Test Kit on nasal swabs, 110 positive specimens and 220 negative specimens confirmed by a CE marked RT-PCR test kit were evaluated.

SARS-CoV-2 Antigen Test Performance against the RT-PCR Method with Nasal Swab Specimen

SARS-CoV-2 Antigen Test Kit	Comparator Method (RT-PCR)	Positive	Negative	Total

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